




 **FirstCare**
A Persys Medical Company



Email: info@firstcareproducts.com | Web: www.firstcareproducts.com

MULTI Bandage Data Sheet

Making the Essential Simple

THE MULTI BANDAGES

The Multi bandage line provides large bandages in compact packaging for treating Abdominal and exposed organ wounds or amputation injuries. Multi bandages are multifunctional with an optional moisture seal for treating exposed organs. Multi bandages are compact, light-weight and vacuum-sealed, suitable for first aid combat emergency kits. The Multi Bandage features a large sterile dressing pad, moisture seal, an elastic cloth wrapping leader and an easy to use closer bar – all in a compact, waterproof packaging.

MULTI BANDAGE FEATURES

- Multi-Functional bandage for large wounds
- Large sterile pad - effective dressing for extensive wounds, burns or amputations
- Moisture Seal for covering exposed organs, this seal can be easily removed if required
- Compact, Light-Weight, Waterproof, Vacuum-sealed packaging
- Sterile Bandage with an 8 year shelf life

Learn more at: www.firstcareproducts.com

ABOUT FIRST CARE PRODUCTS

First Care Products develops and manufactures innovative first-response medical products to improve field hemorrhage control and to save lives. The First Care Emergency Bandage has been adopted by military and civilian organizations worldwide.



PRODUCT SPECIFICATIONS

The Multi bandage line provides life saving solutions for treating large wounds.



Type. FCP-19

Description. 4" Abdominal / Large Wound Amputation Dressing
Elastic Bandage Size.

370cm+- 20cm length (when stretched) x 10cm width

Pad Size. 30cm width x 30cm length

Bag Size. 12.5cm width x 24cm length

Moisture Seal. 40cm width x 40cm length

Part No. 101929000



Type. FCP-29

Description. 6" Abdominal / Large Wound Amputation Dressing
Elastic Bandage Size.

370cm+- 20cm length (when stretched) x 15cm width

Pad Size. 30cm width x 30cm length

Bag Size. 12.5cm width x 24cm length

Moisture Seal. 40cm width x 40cm length

Part No. 102929000



Type. 9T

Description. 8" Abdominal / Large Wound Amputation Dressing
Elastic Bandage Size.

500cm+- 20cm length (when stretched) x 20cm width

Pad Size. 30cm width x 30cm length

Bag Size. 16.5cm width x 30cm length

Moisture Seal. 40cm width x 40cm length

Part No. 105929000

